



PRODUCT CODE: 251301

Eosin Yellowish hydroalcoholic solution 1% (CE-IVD) for clinical diagnostics

Solution for Hematoxylin-Eosin staining, in human and gynaecological samples.

,In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3204 12 00 90

PHYSICAL DATA: Clear , red liquid with green fluorescence if it is observed on a dark background, Miscible with alcohol and water • D 20/4 0,992 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l-0,997 kg 1kg-1,003 l

COMPOSITION:

Eosin Yellowish	1 g
Ethanol absolute	10 ml
Water	90 ml

SPECIFICATIONS:

Identity :	
Identity	passes test

Maximum limit of impurities

Suitability for hematoxylin - eosin stain	passes test
---	-------------

