



**PRODUCT CODE: 255298**

**Carazzi's Hematoxylin solution (CE-IVD) for clinical diagnostics**

Solution for Hematoxylin-Eosin staining, in human and gynaecological samples., In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

---

TARIC 3822 90 00 00

**PHYSICAL DATA:** liquid, Clear, Violet, Miscible with water • D 25/4 1,08 •

**BIBLIOGRAPHY:** Eduardo Suarez Peregrin., "Manual Técnico de Análisis Clínicos", 9 Ed.,Prieto, 388 (1972) ; J.D. Bancroft, A. Stevens, "Theory and Practice ofHistological Techniques", Edinburgh, Scotland •

---

**HAZARDOUS:**

---

**WEIGHT/VOLUME INFORMATION:** 1l-1,080 kg 1kg-0,926 l

**OBSERVATIONS:** Filter any sediment present in the solution. •

---

**COMPOSITION:**

Hematoxylin	0,1 g
Aluminium Potassium Sulphate 12-hydrate	5 g
Sodium Iodate	0,02 g
Glycerol	20 ml
Water (s.q.m)	100 ml

**SPECIFICATIONS:**

Identity :  
Identity passes test

**Maximum limit of impurities**

Suitability for hematoxylin - eosin stain passes test

