



PRODUCT CODE: 255811

Mounting Medium for substitutes of xylene for clinical diagnostics
for microscopy

TARIC 3822 90 00 00

Acrylic modified resin dissolved in isoparaffin H .

PHYSICAL DATA: liquid, Clear, Colourless, Inmiscible with water. • D 20/4 0,82 • Flash P.:63 °C • Vap. press. (20 °C) 0,7 hPa • Expl. limit 0,6 %(V) 7,0 %(V) •

INSTRUCTIONS FOR USE: After processing the slide in the last Isoparaffin tray, place carefully three separated drops of Mounting Medium for substitutes of xylene on the slide.

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat >10.000 mg/kg •



H: H317 • H304 • H413 •

P: P261 • P301+P310 • P321 • P405 • P362+P364 • P331 • P501 •

WEIGHT/VOLUME INFORMATION: 1l~0,820 kg 1kg~1,22 l

SPECIFICATIONS:

Density at 25°C

0,815-0,825

