



Karl Fischer water standards

Water Standards

They are used to calculate the titre of **AQUAMETRIC®** Karl Fischer reagents for volumetric titration and to check the reliability of water determination using the coulometric method.

They are NIST traceable.

The new package

The glass ampoule has been replaced by a glass container closed with a septum which allows the extraction of the contents with a syringe.

Main advantages

- **Safety has been improved**
 - » The user does not come into contact with the product.
 - » Less risk of injury in the laboratory: opening the container is safer compared to glass ampoules.
 - » This new container is not as fragile as the glass ampoules.
- **Results more reliable**
 - » During the analysis the exposure of the product to atmospheric moisture is minimized.
- **No changes in the product stability:** 2 years shelf-life.



Product code	Product name	Volumetric	Coulometric	Pack size
395459.2522	Karl Fischer Water Standard 1.0 mg/g (Reag. Ph. Eur.)		•	10 x 10 mL
395458.2522	Karl Fischer Water Standard 10 mg/g (Reag. Ph. Eur.)	•		10 x 10 mL



Distributors

For ordering our high-quality products please contact your local distributor. **ITW Reagents** has distributors in the countries as follows:

Algeria	Egypt	Lithuania	Slovenia
Argentina	Finland	Mexico	South Korea
Armenia	France	Montenegro	Spain
Australia	Germany	Morocco	Sweden
Austria	Greece	Netherlands	Switzerland
Bangladesh	Hong Kong-China	New Zealand	Taiwan
Belarus	Hungary	Norway	Thailand
Belgium	Iceland	Pakistan	Tunisia
Bosnia and Herzegovina	India	Paraguay	Turkey
Brazil	Indonesia	Peru	Ukraine
Bulgaria	Ireland	Poland	United Kingdom
Chile	Israel	Portugal	Uruguay
China	Italy	Romania	United States
Colombia	Ivory Coast	Russia	Vietnam
Croatia	Japan	Saudi Arabia	
Czech Republic	Jordan	Senegal	
Denmark	Kazakhstan	Serbia	
Ecuador	Kuwait	Slovakia	

For further contact details of our distributors please visit

itwreagents.com

or contact our Customer Service at

customerservice.de@itwreagents.com

customerservice.es@itwreagents.com

customerservice.it@itwreagents.com

Customers - Webshop

Authorized customers from Spain, Portugal, Austria, France and Germany can order via our webshop on

itwreagents.com.

IP-042EN

AppliChem GmbH
Ottoweg 4
DE-64291 Darmstadt
Germany
Phone +49 6151 9357 0
Fax +49 6151 9357 11
info.de@itwreagents.com

Nova Chimica Srl
Via G. Galilei, 47
I-20092 Cinisello Balsamo
(Milano) Italy
Phone +39 02 66045392
Fax +39 02 66045394
info.it@itwreagents.com

PanReac Química SLU
C/ Garraf 2, Polígono Pla de la Bruguera
E-08211 Castellar del Vallès
(Barcelona) Spain
Phone +34 937 489 400
Fax +34 937 489 401
info.es@itwreagents.com



www.itwreagents.com