


PRODUCT CODE: 132440

Sodium tetra-Phenylborate (Reag. USP, Ph. Eur.) for analysis, ACS
 (KALIGNOST®) (® Registered trade-mark of Heyl Co. Berlin)

C₂₄H₂₀BNa

C₂₄H₂₀BNa



M.= 342,23

CAS [143-66-8]

EINECS 205-605-5

TARIC 2931 90 00 90

SYNONYMS: Kalignost (Registered trade-mark of Heyl Co.Berlin), Sodium Tetraphenylborate, Tetraphenylboron Sodium

PHYSICAL DATA: crystals, White, Soluble in water. D 0,5 • pH(50 g/l)8 •

BIBLIOGRAPHY: Merck Index **12**, 8.839 13, 8.765 Römp **8** , **2747** • Beilstein **16** , **IV**, **1624** • BRN 3599783 • ACS **XI** •

HAZARDOUS:

SPECIFICATIONS:

Minimum assay	99,5%
Identity :	
Identity	passes test
Maximum limit of impurities	
Appearance of solution	passes test

Loss on drying at 105°C
Sensitivity to K

0,5%
passes test