


**PRODUCT CODE: 172417**

**o-Tolidine solution 0.1%**

Chlorine reagent

C<sub>14</sub>H<sub>16</sub>N<sub>2</sub>

C<sub>14</sub>H<sub>16</sub>N<sub>2</sub>



M.= 212,29

CAS [119-93-7]

TARIC 2921 59 50 00

**SYNONYMS:** 3,3'-Dimethyl-(1,1'-Biphenyl)- 4,4'-Diamine, 3,3'-Dimethylbenzidine

**PHYSICAL DATA:** liquid, Clear, Colourless, Miscible with water • D 20/4 1,021 • pH1 •

**BIBLIOGRAPHY:**

**HAZARDOUS:** LC50 rat 3124ppm (V) / 1h • VLA-EC (HCl) 10 ppm15 mg/m<sup>3</sup> VLA-ED (HCl) 5 ppm7,6 mg/m<sup>3</sup>



H: H350 • H290 •

P: P234 • P280 • P308+P313 • P406 • P405 • P501 •

**WEIGHT/VOLUME INFORMATION:** 1l-1,021 kg 1kg-0,979 l

**COMPOSITION:**

o-Tolidine	0,1 g
Hydrochloric Acid 35%	10,5 ml
Water (s.q.m)	100 ml

