



PRODUCT CODE: 175305

Sodium Chloride solution ASTM B117, ISO 9227

for saline spray (fog)

NaCl

NaCl

M.= 58,44

CAS [7647-14-5]

EINECS 231-598-3

TARIC 2501 00 99 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,035 •

BIBLIOGRAPHY: UNE-EN ISO 9227:2017 Corrosion testing in artificial atmospheres. Sal mist testing. • ASTM Standard B117-19 Standard Practice for Operating Salt Spray (Fog) Apparatus, UNE-EN 60068-2-11 Environmental testing. Part 2: Tests. Test Ka: Salt mist •

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,035 kg 1kg~0,966 l

SPECIFICATIONS:

Assay (Arg.)	4,5-5,5 %
pH at 23+-3 / 35 °C	6,5-7,2
Density at 25/4	1,029 - 1,036

Maximum limit of impurities

Anti-caking Agents	0,0%
Sodium Iodide (NaI)	0,1 %
Halides (Bromide, Fluoride, Iodide)	0,1%
Total impurities	0,2%

Cu	0,00003 %
Ni	0,001 %
Pb	0,0025%

This product has been manufactured with Sodium Chloride which passes the ASTM B117-19 and ISO 9227-17 specifications and Water Type IV which passes ASTM D1193-06 specifications.

Ed.: 9 . Vig.: 21.02.2020 .

Prod.: 175305