



PRODUCT CODE: 211630

Silicone paste A technical grade
for greasing at high temperatures

TARIC 3910 00 00 90

Mixture of polyorganosiloxanes and inert loads .

PHYSICAL DATA: Soft yellowish-white mass, Insoluble in water Thermal decom. > 200 °C • D (20 °C) 1,1 • Flash P.:>400 °C •

BIBLIOGRAPHY:

HAZARDOUS: LD50 po rat > 5.000 mg/kg LD50 pc rat > 2.000 mg/kg TWA (SiO₂) 10mg/m³/8h

SPECIFICATIONS:

Identity :
Identity

IR passes test