



PRODUCT CODE: 211631

Silicone paste B technical grade

for greasing at pressure and vacuum. For use at pressures down to 10⁻⁶ mm Hg.

TARIC 3910 00 00 90

Mixture of polyoxydimethylsiloxanes and inert loads .

PHYSICAL DATA: Sticky mass, Greyish-white, Insoluble in water Thermal decom. > 200 °C • D (25 °C) 1,01 •
Flash P.: >400 °C •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat >5.000 mg/kg • LD50 po rat > 5.000 mg/kg LD50 pc rat > 2.000 mg/kg

SPECIFICATIONS:

Identity :
Identity

IR passes test