



PRODUCT CODE: 216241

Silicone antifoaming liquid (AQ) technical grade
for distillation of aqueous products

TARIC 3910 00 00 90

10% Active silicone emulsion .

PHYSICAL DATA: liquid, White, D 20/4 1 • B.P.:> 35 °C • pH25 °C3,5 • Flash P.:>100 °C •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,00 kg 1kg~1,00 l

SPECIFICATIONS:

Identity :
Identity

IR passes test