



**PRODUCT CODE: 251172**

**Kühne's Methylene Blue Phenicated solution (CE-IVD) for clinical diagnostics**  
for microscopy, In Vitro Diagnostic medical device class A in compliance to the  
Regulation (EU) 2017/746.

---

TARIC 2934 99 90 90

**PHYSICAL DATA:** liquid, Dark blue, Miscible with water D 20/4 0,995 •

**BIBLIOGRAPHY:** Manual Técnico de Análisis Clínico, 9 Ed, E. Suárez, 171 (1972) •

---

**HAZARDOUS:**



H: H332 • H319 • H315 • H341 •

P: P264 • P270 • P301+P312 • P302+P352 • P312 • P501 • P322 • P330 • P363 • P308+P313 •

---

**WEIGHT/VOLUME INFORMATION:** 1l-0,995 kg 1kg-1,005 l

---

**COMPOSITION:**

Methylene Blue	9 g
Ethanol absolute	90 ml
Phenol	26 ml
Water	1000 ml

**SPECIFICATIONS:**

Identity :  
Identity passes test

**Maximum limit of impurities**

Staining suitability

passes test

---

Ed.: 4 . Vig.: 19.05.2023 .

Prod.: 251172