



**PRODUCT CODE: 252164**

**Buffer Solution pH 7.2 for clinical diagnostics**  
for hematology, according to Weise

---

TARIC 3822 90 00 00

**PHYSICAL DATA:** liquid, Clear, Colourless, Miscible with water • D 20/4 1,001 • pH7,2 •

**BIBLIOGRAPHY:**

---

**HAZARDOUS:**

---

**WEIGHT/VOLUME INFORMATION:** 1l~1,001 kg 1kg~0,999 l

---

**COMPOSITION:**

Potassium di-Hydrogen Phosphate	40 mg
di-Sodium Hydrogen Phosphate 12-hydrate	151 mg
Water (s.q.m)	100 ml

**SPECIFICATIONS:**

Identity :	
Identity	passes test
pH	7,0 - 7,4