



PRODUCT CODE: 253998

Blue for fast staining (Panoptic No. 3)(CE-IVD) for clinical diagnostics

blood smear staining or medullary smear staining, In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Dark blue, Miscible with water D 20/4 1,007 • pH6,8 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l-1,007 kg 1kg-0,993 l

COMPOSITION:

Azur B	2 g
Buffer solution pH 7 s.q.m	1000 ml

SPECIFICATIONS:

Identity :	
Identity	passes test

Maximum limit of impurities

Suitability for hematology stain	passes test
----------------------------------	-------------