



PRODUCT CODE: 253999

Eosin for fast staining (Panoptic No. 2)(CE-IVD) for clinical diagnostics

blood smear staining or medullary smear staining, In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Orange-red, Miscible with water D 20/4 1,006 • pH6,8 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l-1,006 kg 1kg-0,994 l

COMPOSITION:

Eosin Yellowish	0,8 g
Buffer solution pH 7 s.q.m	1000 ml

SPECIFICATIONS:

Identity :	
Identity	UV passes test
pH	7,00 ± 0,05

Maximum limit of impurities

ABS	
at λ 517 nm; dil. 1:100	0,948 - 1,048
Staining suitability	passes test

