



PRODUCT CODE: 256993

Paraffin M.P. 55-58°C plasticized + DMSO pellets(CE-IVD) for clinical diagnostics

Embedding medium for preparation of fixed and dehydrated human tissue samples., In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

CAS [8002-74-2]

EINECS 232-315-6

TARIC 2712 20 90 00

PHYSICAL DATA: pellets, White, Insoluble in water M.P.: 55 - 58 °C • B.P.:> 350 °C • Ign. T.:245 °C •

BIBLIOGRAPHY: Römp 8 , 2991 • Royal Decree I1418 •

HAZARDOUS: VLA-ED 2 mg/m3

SPECIFICATIONS:

Identity :

Identity

Melting range

IR passes test

55 - 58 °C

Maximum limit of impurities

Acidity or alkalinity

Insoluble matter in C8H10

passes test

passes test