



PRODUCT CODE: 257462

Histofix ® Preservative ready to use (pink)(CE-IVD) for clinical diagnostics

In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3822 90 00 00

PHYSICAL DATA: Clear liquid, Pink, Miscible with water D 20/4 1,019 • pH6,8 - 7,2 • It is decolorized with exposure to light

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 100 mg/kg • LC50 rat 203mg/m³ • LD50 skn rbt 270 mg/kg • VLA-EC (Formaldehyde) 0,37 mg/m³ VLA-EC (Formaldehyde) 0,3 ppm



H: H350 • H317 • H341 • H302 •

P: P201 • P202 • P261 • P272 • P280 • P281 • P302+P352 • P308+P313 • P321 • P333+P313 • P363 • P405 • P501 •

WEIGHT/VOLUME INFORMATION: 1l-1,019 kg 1kg-0,981 l

OBSERVATIONS: Storage away from direct light. •

SPECIFICATIONS:

Assay (Iodom.)	3,7 - 4,0 %
pH	6,8 - 7,2
UV Spectrum (1cm cell; Ref.: water) ABS λ máx. (525 - 527 nm)	>= 0.9

Methanol (w/v)

1 - 1,5 %

Ed.: 4 . Vig.: 22.12.2023 .

Prod.: 257462