



PRODUCT CODE: 273526

TISAB III - Solution for detection of fluoride

for the fluoride determination by selective electrodes. pH 5,25 ±0,25

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Miscible with water • D 20/4 1,053 • pH5,0 - 5,5 •

INSTRUCTIONS FOR USE: Add TISAB III to samples and standards in a ratio 1:9.

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l-1,053 kg 1kg-0,95 l

COMPOSITION:

1,2-Diaminociclohexane N,N,N,N-tetraacetic Acid 1-hydrate	18 g
Ammonium Chloride	96,65 g
Ammonium Acetate	163,4 g
Cresol Red	0,1 g
Water (s.q.m)	1 l

Note: pH 5.25±0.25