

Specification

**Buffer Solution pH 12.454**

**277089**

<b>Physical Description:</b>	Liquid
<b>Product Code:</b>	277089
<b>Product Name:</b>	Buffer Solution pH 12.454
<b>Short Description:</b>	usage: These buffer solutions are useful for the calibration of high precision laboratory pH meters
<b>Headline Comment:</b>	* standardized, according to DIN19266
<b>Specifications:</b>	<p>pH (10°C): 12.810 ± 0.01</p> <p>pH (20°C): 12.627 ± 0.01</p> <p>pH (25°C): 12.454 ± 0.01</p> <p>pH (30°C): 12.289 ± 0.01</p>
<b>WGK:</b>	nwg
<b>Storage:</b>	RT
<b>CS:</b>	38221900

**AppliChem GmbH**

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • [info.de@itwreagents.com](mailto:info.de@itwreagents.com) • [www.itwreagents.com](http://www.itwreagents.com)  
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340