

Specification

Buffer Solution pH 6.865

277124

Physical Description:	Liquid
Product Code:	277124
Product Name:	Buffer Solution pH 6.865
Short Description:	usage: These buffer solutions are useful for the calibration of high precision laboratory pH meters
Headline Comment:	* standardized, according to DIN19266
Specifications:	<p>pH (10°C): 6.900 ± 0.01</p> <p>pH (20°C): 6.881 ± 0.01</p> <p>pH (25°C): 6.865 ± 0.01</p> <p>pH (30°C): 6.853 ± 0.01</p>
WGK:	nwg
Storage:	RT
CS:	38221900

AppliChem GmbH