



PRODUCT CODE: 281384

Acidimetric Liquor titrated for volumetric analysis

for determination of acidity in milk. 1 ml corresponds to 0,01g of lactic acid (=1° Dornic)

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,006 • pH13 •

INSTRUCTIONS FOR USE: Stir 10 ml of milk with 10 drops of phenolphthalein 1%. Titrate using a graduated burette with Acidimetric Liquor stirring constantly until a persistent pink colour appears. The ml consumed in the titration correspond to °Dornic.

BIBLIOGRAPHY: Métodos Oficiales de Análisis IV , 60 and 61 (1993) •

HAZARDOUS: LD L0 oral rbt 500 mg/kg • LD50 ipr mus 40 mg/kg

WEIGHT/VOLUME INFORMATION: 1l~1,006 kg 1kg~0,994 l

SPECIFICATIONS:

Titer at 20°C	1,000±0,005
Sodium Hydroxide	0,111 mol/l