


PRODUCT CODE: 281440

Eriochrome Black T solution 1% for volumetric analysis
 for complexometry

$C_{20}H_{12}N_3NaO_7S$

$C_{20}H_{12}N_3NaO_7S$



M.= 461,38

TARIC 3822 90 00 00

SYNONYMS: 3-Hydroxy-4-[(1-Hydroxy-2-Naphtalenyl)azo]-7-Nitro-1-Naphtalenesulphonic Acid Sodium Salt, Mordant Black 11

PHYSICAL DATA: liquid, viscous, Dark blue, Miscible with alcohol and water • D 20/4 1,05 •

BIBLIOGRAPHY:

HAZARDOUS:



H: H225 • H314 • H318 •

P: P310 • P303+P361+P353 • P305+P351+P338 • P501 • P280 • P321 • P405 •

TRANSPORT REGULATIONS: UN: 1993 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 364 • CAO: 353 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~1,050 kg 1kg~0,952 l

COMPOSITION:

Eriochrome Black T	1 g
Triethanolamine	75 ml
Ethanol absolute	25 ml

