




PRODUCT CODE: 363622

Tetrabutylammonium Hydrogen Sulfate for HPLC
for ion pair chromatography

$C_{16}H_{37}NO_4S$

$C_{16}H_{37}NO_4S$



M.= 339,54 CAS [32503-27-8] EINECS 251-068-5

TARIC 2923 90 00 90

SYNONYMS: Tetrabutylammonium bisulphate, tetra-n-Butylammonium Hydrogen Sulphate

PHYSICAL DATA: Crystalline powder or Crystals, White to slightly yellowish, Soluble in water. Hygroscopic · M.P.: 169 - 171 C ·

BIBLIOGRAPHY: Safety 2 , 3237 D · Fieser 6565 7354 8461 9356 1041 11403 12380 1394 1515 · Lett. 29 , 4567 (1988) ·

HAZARDOUS:



H: H315 · H319 ·

P: P264 · P280 · P305+P351+P338 · P332+P313 · P337+P313 · P501 ·

SPECIFICATIONS:

Minimum assay (Acidim.)

99%

Identity :

Identity

UV Spectrum (1cm cell; Ref.: water):

Transmittance at 200 nm

Transmittance at 220 nm

IR passes test

? 70 %

? 89 %

Ed.: 4 . Vig.: 06.03.2015 .

Prod.: 363622