

Specification

Tetrabutylammonium di-Hydrogen Phosphate for HPLC

367038

Physical Description:	Solid
Product Code:	367038
Product Name:	Tetrabutylammonium di-Hydrogen Phosphate for HPLC
Specifications:	<p>Assay (titr., calc. on dried subst.): min. 99 %</p> <p>pH (1 %; H₂O; 20°C): 3.0 - 5.0</p> <p>Loss on drying (2 h; 110°C): max. 6 %</p> <p>A (1 cm/0.005 M in H₂O)</p> <p>210 nm: max. 0.15</p> <p>230 nm: max. 0.05</p> <p>260 nm: max. 0.02</p>
WGK:	1
Storage:	RT
Molecular Formula:	C ₁₆ H ₃₈ NO ₄ P
M:	339.46 g/mol
CAS:	5574-97-0
EINECS:	226-947-1
CS:	29239000

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340