



PRODUCT CODE: 394658

Conductivity Standard 12.88 mS/cm (25°C)

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Miscible with water • D 1,002 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l-1,002 kg 1kg-0,998 l

Potassium Chloride 0,100 mol/l

SPECIFICATIONS:

CONDUCTIVITY TEMPERATURE: T (°C)	VALUES	DEPENDING	ON
20,0			(mS/cm) 11,67
21,0			11,91
22,0			12,15
23,0			12,39
24,0			12,64
25,0			12,88
26,0			13,13
27,0			13,37
28,0			13,62
29,0			13,87
30,0			14,12

Uncertainty

See certificate

Conductivity at 25°C

12,82 - 12,94 mS/cm