



PRODUCT CODE: 395443

Redox Standard 220 mV (25°C)

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, yellow or greenish yellow, Miscible with water • D 20/4 1,004 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,004 kg 1kg~0,996 l

COMPOSITION:

Potassium Hexacyanoferrate(II) 3-hydrate	4,3 g
Potassium Hexacyanoferrate(III)	3,3 g
Buffer Solution pH 7.00	200 ml
Water (s.q.m)	1000 ml

SPECIFICATIONS:

REDOX POTENTIAL VALUES DEPENDING ON TEMPERATURE:

T (°C)	mV
10	245
15	236
20	228
25	220
30	212
35	204
40	195
50	178
60	160
70	142

Uncertainty

See certificate

Potencial at 25°C

215 -225 mV