



**PRODUCT CODE: 395508**

**Colour Standard Pt-Co, 500 APHA**  
for colour determination

---

TARIC 2806 10 00 00

**SYNONYMS:** Hazen Colour

**PHYSICAL DATA:** liquid, Golden yellow, Miscible with water • D (20 °C) 1,02 •

**BIBLIOGRAPHY:**

**HAZARDOUS:**



H: H350i • H290 •

P: P201 • P202 • P280 • P308+P313 • P405 • P501 •

**WEIGHT/VOLUME INFORMATION:** 1l~1,020 kg      1kg~0,98 l

**OBSERVATIONS:** Product controlled as a drug precursor. •

**COMPOSITION:**

Potassium Chloroplatinate	1,245 g
Cobalt (II) chloride 6-H <sub>2</sub> O	1,000 g
Water (s.q.m)	1000 ml

**SPECIFICATIONS:**

Uncertainty of the Absorbance	See certificate
-------------------------------	-----------------

ABS

at  $\lambda$  510 nm

at  $\lambda$  430 nm

at  $\lambda$  480 nm

at  $\lambda$  455 nm

0,055 - 0,065

0,110 - 0,120

0,105 - 0,120

0,130 - 0,145

Ed.: 3 . Vig.: 10.07.2014 .

Prod.: 395508