



PRODUCT CODE: 403682

Tryptone (Ingredient) for microbiology

Physical Data: Dehydrated product for culture media of cream colour ·Description: The tryptone is obtained by pancreatic degradation of casein. It is used as a source of nitrogen in culture media for detection of fungi and some bacteria ·

TARIC 3504 00 90 00

PHYSICAL DATA: Soluble in water. Dehydrated product for culture media of cream colour•

BIBLIOGRAPHY:

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place. •

SPECIFICATIONS:

pH sol. 2%	6,5-7,5
Loss on drying at 105°C	6%
Residue on ignition	15 %
Nitrogen, total	>= 10 %