



**PRODUCT CODE: 403687**

**Yeast Extract (Ingredient) for microbiology**

Physical Data: Dehydrated product for culture media of beige colour ·Description: The yeast extract is obtained from autolyzed yeast cells and is rich in vitamin B complex ·

---

TARIC 3504 00 90 00

**PHYSICAL DATA:** Dehydrated product for culture media of beige colour, in water at 20°C D 0,5 • pH(10 g/l sol.)7,0 • Ign. T.:>200 °C •

**BIBLIOGRAPHY:**

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**SPECIFICATIONS:**

pH sol. 2%	6,0-7,2
Dry matter	>= 94 %
Nitrogen, total	>=10 %