



PRODUCT CODE: 403695

Peptone, Bacteriological (Ingredient) for microbiology

Physical Data: Ingredient base for culture media of yellowish colour Description: The Bacteriological Peptone is used as a source of nitrogen for cultivating less fastidious microorganisms. It's water soluble at the usual working concentrations of the culture media Storage in a cool and dry place

EINECS 293-428-4
TARIC 3504 00 90 00

PHYSICAL DATA: Soluble in water.

BIBLIOGRAPHY:

HAZARDOUS:

SPECIFICATIONS:

| | |
|------------------------|---------|
| pH sol. 2% | 6,5-7,5 |
| Loss on drying at 105C | 6% |
| Residue on ignition | 15 % |
| Nitrogen, total | >=10 % |