



PRODUCT CODE: 403898

Casein Peptone (Ingredient) for microbiology

Physical Data: Fine powder of yellow colour ·Description: The Casein Peptone is a pancreatic digestion of the casein with high tryptophan content. It dissolves easily and produces very transparent solutions ·

TARIC 3504 00 90 00

PHYSICAL DATA: Fine powder, Yellow, in water at 20°C pH(20 g/l sol.)6 - 7 ·

BIBLIOGRAPHY:

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

SPECIFICATIONS:

pH sol. 2%	6,5-7,5
Loss on drying at 105°C	6%
Residue on ignition	15 %
Nitrogen, total	>=10 %