



PRODUCT CODE: 413745

Violet Red Bile Glucose Agar (VRBG) (Ph. Eur.) (ISO21528) (Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of reddish beige colour · Preparation: Prepare a suspension of 41,5 g in 1 l of distilled water. Heat to boiling point while constantly stirring and boil for about 1 minute. Let it cool to 45°C and distribute. Alternatively sterilize at 118°C for 15 minutes. AVOID OVERHEATING OR REMELT THE MEDIUM ·

TARIC 3821 00 00 00

SYNONYMS: Crystal Violet-Neutral Red-Bile-Glucose Agar according to Mossel, VRBD Agar, VRBG Agar•

PHYSICAL DATA:

BIBLIOGRAPHY: ISO 21528 • USP 38 • Ph. Eur. 8.0 • J. Appl. Bact., 26, 444-452 (1963), J. Bact., 84, 381 (1962)

•

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Bile Salts Mixture	1,5
Crystal Violet	0,002
Neutral Red	0,03
D(+)-Glucose	10,0
Yeast Extract	3,0
Pancreatic Digest of Gelatine	7,0
Sodium Chloride	5,0
Agar	15,0
pH: 7.4±0.2	

Ed.: 4 . Vig.: 27.06.2011 .

Prod.: 413745