



PRODUCT CODE: 413777

Brain Heart Infusion (BHI) (Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of beigecolour · Preparation: Make a suspension of 37 g in 1 l of distilled water. Heat stirring constantly to boiling point and boil for 1 minute. Distribute and sterilize at 121°C for 15 minutes ·

TARIC 3821 00 00 00

SYNONYMS: Brain Heart Broth

PHYSICAL DATA:

BIBLIOGRAPHY: J. Bact., 62 , 613 (1951) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place. •

Composition (g/l):

Calf Brain Infusion	7,5
Meat Heart Infusion	10,0
D(+)-Glucose	2,0
Gelatin Peptone	10,0
Sodium Chloride	5,0
di-Sodium Hydrogen Phosphate	2,5
pH: 7.4±0.2	

