



**PRODUCT CODE: 413780**

**MacConkey Broth (Ph. Eur.) (Dehydrated Culture Media) for microbiology**

Physical Data: Dehydrated culture media of beige colour ·Preparation: Dissolve 35 g in 1 l of distilled water. If the sample to test is bigger than 1 ml, prepare medium double strength. Distribute into test tubes containing Durham bell and sterilize at 121°C for 15 minutes ·

---

TARIC 3821 00 00 00

**SYNONYMS:** Medium G·

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** USP 38 · Ph. Eur. 8.0 · J. Hyg., 8, 322-334 (1908) ·

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

Ox Bile	5,0
Lactose	10,0
Gelatin Peptone	20,0
Bromocresol Purple	0,01
pH: 7.3±0.2	