



PRODUCT CODE: 413781

Malt Extract Agar (Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of whiteto light beige colour · Preparation: Make a suspension of 33,6 g in 1 l of distilled water. Heat stirring constantly to boiling point and boil for 1 minute. Sterilize at 118°C for 10 minutes. If the medium is overheated, the agar loses its capacity to solidify ·

TARIC 3821 00 00 00

PHYSICAL DATA:

BIBLIOGRAPHY: Thom & Church. The Aspergilli. (1926) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Malt Extract	12,75
Dextrin	2,75
Glycerol	2,35
Gelatin Peptone	0,78
Agar	15,0
pH: 4.7±0.2	