



PRODUCT CODE: 413785

MRS Broth (Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of creamcolour · Preparation: Make a suspension of 52 g in 1 l of distilled water. Heat stirring until dissolved. Distribute and sterilize at 121°C for 15 minutes ·

TARIC 3821 00 00 00

SYNONYMS: De Man, Rogose, Sharpe, Broth•

PHYSICAL DATA:

BIBLIOGRAPHY: J. Appl. Bact., 22, 329-340 (1959), J. Appl. Bact., 23, 130-135 (1960) •

HAZARDOUS:

OBSERVATIONS: Storage between +2 and +8°C, in a dry place, away from direct light.

Composition (g/l):

di-Ammonium Hydrogen Citrate	2,0
Meat Extract	8,0
Yeast Extract	4,0
D(+)-Glucose	20,0 g
Magnesium Sulfate	0,2
Manganese(II) Sulphate	0,05
Bacteriological Peptone	10,0
di-Potassium Hydrogen Phosphate	2,0
Sodium Acetate	5,0
Tween 80	1,0
pH: 6.2±0.2	

Ed.: 3 . Vig.: 15.01.2011 .

Prod.: 413785