



PRODUCT CODE: 413788

Mueller-Hinton Broth (Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of beigecolour · Preparation: Make a suspension of 21 g in 1l of distilled water.Heat stirring constantly to boiling point and boil for 1 minute. Distribute and sterilize at 121°C for 15 minutes. Do not overheat ·

TARIC 3821 00 00 00

PHYSICAL DATA:

BIBLIOGRAPHY: Amer. J. Clin. Pathol., 45: 493-496 (1966), J. Clin. Microbiol., 22: 369-374 (1985), Standardization of Methods for Conducting Microbic Sensitivity Test. WHO. (1961) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Starch	1,5
Meat Infusion	2,0
Casein Peptone Hydrolysate	17,5
pH: 7.4±0.2	