



PRODUCT CODE: 413792

Nutrient Agar (ISO 6579, ISO 10273, ISO 19250)(Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of light beige colour · Preparation: Make a suspension of 23 g in 1 l of distilled water. Heat stirring constantly to boiling point and until completely dissolved. Boil for 1 minute. Sterilize at 121°C for 15 minutes ·

TARIC 3821 00 00 00

PHYSICAL DATA:

BIBLIOGRAPHY: Standard Methods for the Examination of Water and Wastewater. APHA. (1980). (ISO 6579, ISO 10273, ISO 19250)., Atlas R.M. & Parks L.C., Handbook of Microbiological Media, 667 (1993) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Meat Extract	3,0
Gelatin Peptone	5,0
Agar	15,0
pH: 6.8±0.2	