



**PRODUCT CODE: 413802**

**Sabouraud Glucose Agar (Ph. Eur.) (Dehydrated Culture Media) for microbiology**

Physical Data: Dehydrated culture media of lightbeige colour · Preparation: Prepare a suspension of 65 g in 1 l of distilled water. Heat stirring constantly to boiling point and boil for 1 minute. Sterilize at 118-121°C for 15 minutes ·

---

TARIC 3821 00 00 00

**SYNONYMS:** Sabouraud Agar, Sabouraud-4 %-Glucose Agar·

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** USP 38 · Ph. Eur. 8.0 ·

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

D(+)-Glucose	40,0
Agar	15,0
Mixture of Peptic Digest of Animal Tissue and Pancreatic Digest of Casein (1:1)	10
pH: 5.6±0.2	