

**PRODUCT CODE: 413805****Salmonella Shigella Agar (Dehydrated Culture Media) for microbiology**

Physical Data: Dehydrated culture media of lightyellow to pink colour · Preparation: Make a suspension of 60 g in 1 l of distilled water. Heat stirring constantly to boiling point and boil for 1 minute. DO NOT STERILIZE IN AUTOCLAVE. Allow to cool to 45-50°C and distribute into sterile Petri dishes, using 20 ml per dish. Leave the medium to solidify ·

---

TARIC 3821 00 00 00

**SYNONYMS:** SS Agar•**PHYSICAL DATA:****BIBLIOGRAPHY:** Atlas R.M. & Parks L.C., Handbook of Microbiological Media , 788 (1993), Compendium of Methods for the Microbiological Examination of Foods 2nd ed. APHA. (1984) •**HAZARDOUS:****OBSERVATIONS:** Storage in a cool and dry place.**Composition (g/l):**

Meat Extract	5,0
Iron(III) Citrate	1,0
Lactose	10,0
Peptones	5,0
Neutral Red	0,025
Bile Salts	8,5
tri-Sodium Citrate	8,5
Sodium Thiosulphate	8,5
Brilliant Green	0,00033
Agar	13,5
pH: 7.0±0.2	

Ed.: 1 . Vig.: 15.01.2011 .

Prod.: 413805