



**PRODUCT CODE: 414723**

**Egg Yolk Tellurite Emulsion (Supplement) for microbiology**

Physical Data: Additive, stabilized for bacteriological use ·Preparation: Add, aseptically, 50 ml per litre of melting culture media, sterile and cooled to 55-60°C. It's an additive for: Tryptone Soy Agar (TSA) (Cod. 413819), Nutrient Agar (Cod. 413792) and Baird-Parker Agar Base (Cod. 413744) ·

---

TARIC 3821 00 00 00

**PHYSICAL DATA:**

**BIBLIOGRAPHY:**

**HAZARDOUS:**

**OBSERVATIONS:** Storage between +2 and +8°C, in a dry place, away from direct light.