



PRODUCT CODE: 414724

Potassium Tellurite solution 3.5% (Supplement) for microbiology

Physical Data: Additive, stabilized for bacteriological use · Preparation: Add, aseptically, the Potassium Tellurite solution to the melting and sterile culture media below 45°C.

TARIC 3821 00 00 00

PHYSICAL DATA: Miscible with water

BIBLIOGRAPHY:

HAZARDOUS:

OBSERVATIONS: Storage away from direct light. • Storage in a cool and dry place. •

COMPOSITION:

Potassium Tellurite	35,0 g
Distilled water (s.q.m.)	1000 ml