



PRODUCT CODE: 414753

Luria Broth Base (Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of light beige colour · Preparation: Make a suspension of 25 g in 1 l of distilled water. Heat stirring until completely dissolved. Distribute and sterilize at 121°C for 15 minutes ·

TARIC 3821 00 00 00

SYNONYMS: ATCC Medium 1065, ATCC Medium 1082, LB Broth, Miller, LB Medium •

PHYSICAL DATA:

BIBLIOGRAPHY: Atlas R.M. & Parks L.C., Handbook of Microbiological Media, 492 (1993) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Casein Peptone	10,0
Yeast Extract	5,0
Sodium Chloride	10,0
pH: 7.0±0.2	