

**PRODUCT CODE: 415382****Lethen Broth (modified) (Dehydrated Culture Media) for microbiology**

Physical Data: Dehydrated culture media of tostage beige colour · Preparation: Suspend 43,8 g of the medium in one liter of distilled water. Mix well and dissolve by heating with frequent agitation. Boil for one minute until complete dissolution. Do not overheat. Dispense into appropate containers and sterilize at 121°C for 15 minutes. Storage prepared media between +2 and +8°C. The colour of prepared medium is amber, slightly opalescence. The dehydrated medium should be homogeneous, free-flowing and beige colour. If there are any physical changes, discard the medium ·

TARIC 3821 00 00 00

PHYSICAL DATA:

BIBLIOGRAPHY: Bacteriological Analytical Manual. Chapter 23: Microbiological Methodsfor Cosmetics, Lethen Agar (modified), Lethen Broth (modified). 8thEd. FDA. (1995) •

HAZARDOUS:**Composition (g/l):**

Meat Extract	5,0
Yeast Extract	2,0
Lecithin	0,70
Casein Peptone	5,0
Meat Peptone	20,0 g
Dextrose	1,0
Sodium Chloride	5,0
Sodium Bisulfite	0,10
Polysorbate 80	5,0
pH: 7.2±0.2	

Ed.: 3 . Vig.: 09.12.2014 .

Prod.: 415382