



**PRODUCT CODE: 416106**

**Tryptone Yeast Extract Agar (ISO 6222:1999) (Dehydrated Culture Media) for microbiology**

Physical Data: Dehydrated culture media of beige colour · Preparation: Make a suspension of 24 g in 1 l of distilled water. Heat stirring constantly to boiling point and until completely dissolved. Boil for 1 minute. Sterilize at 121°C for 15 minutes. Cool to 45-50°C and distribute into sterile Petri dishes ·

---

TARIC 3821 00 00 00

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** ISO 6222:1999 Water quality. Enumeration of culturable micro-organisms. Colony count by inoculation in a nutrient agar culture medium ·

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

Yeast Extract	3,0
Tryptone	6,0
Agar	15,0
pH: 7.2±0.2	