



PRODUCT CODE: 433819

Tryptone Soy Agar (TSA) (Ph. Eur.) (Contact Plate) for microbiology

Physical Data: Prepared culture media of beige colour · Procedure: Get the culture medium in contact with the test surface. Incubate at 30-35°C for 24-48 hours. Count the number of colonies, the result as cfu/plate ·

TARIC 3821 00 00 00

SYNONYMS: Medium B, Tryptone Soya Agar, TSA•

PHYSICAL DATA:

BIBLIOGRAPHY: USP 38 • Ph. Eur. 8.0 • J. Clin. Microbiol., 23 , 600-603 (1986), Atlas R.M. & Parks L.C., Handbook of Microbiological Media, 949 (1993) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Soy Peptone	5,0
Casein Peptone	15,0
Sodium Chloride	5,0
Agar	15,0
pH: 7.3±0.2	