



PRODUCT CODE: 445576

TSC Agar (ISO 14189, 7937) (Prepared Plate (@ 55 mm)) for microbiology

Physical Data: Prepared media of beige colour · Procedure: Filter appropriate volume through a sterile 47 mm, 0,45 m m gridded membrane filter, under partial vacuum. Place filter on agar in Petri dish. Incubate in anaerobic conditions for 24 +/-3 hours at 44 +/-1°C. Count the black colonies ·

TARIC 3821 00 00 00

SYNONYMS: FDA M169, Tryptose-Sulphite-Cycloserine Agar

PHYSICAL DATA:

BIBLIOGRAPHY: ISO 6461-2 (1986)ISO 7937 (1997)ISO 14189 • EN26461-2 (1995), FDA Bacteriological Analytical Manual 8 ed. AOAC(1995), ISO 7937, Atlas. R. M. and L. C. Parck (1993). Handbook of Microbiological Media. CRC Press, Inc. London, CeNAN (1982). Técnicas para el análisis Microbiológico de Alimentos y Bebidas. Madrid •

HAZARDOUS:

Composition (g/l):

Yeast Extract	5,0
Iron(III) Citrate	1,0
Soy Peptone	5,0
Sodium Disulphite	1,0
Tryptose	15,0
Cycloserine	0,4
Agar	14,0
pH: 7.6±0.2	

Ed.: 6 . Vig.: 29.11.2018 .

Prod.: 445576