



**PRODUCT CODE: 453792**

**Nutrient Agar (ISO 6579, ISO 10273, ISO 19250)(Prepared Plate (@ 90 mm)) for microbiology**

Physical Data: Prepared media of light beige colour · Procedure: Inoculate the plate with the sample ·

---

TARIC 3821 00 00 00

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** Standard Methods for the Examination of Water and Wastewater. APHA. (1980). (ISO 6579, ISO 10273, ISO 19250)., Atlas R.M. & Parks L.C., Handbook of Microbiological Media, 667 (1993) •

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

Meat Extract	3,0
Meat Peptone	5,0
Agar	12,0
pH: 7.0±0.2	