



PRODUCT CODE: 455095

TSA-Tween-Lecithin-Agar (Ph. Eur.) (Prepared Plate (@ 90 mm)) for microbiology

Physical Data: Prepared plate of beige colour · Procedure: Inoculate the plate with the sample ·

TARIC 3821 00 00 00

SYNONYMS: Tryptone Soy Agar + Tween + Lecithin, TSA-Polysorbate-Lecithin-Agar·

PHYSICAL DATA:

BIBLIOGRAPHY: USP 38 · Ph. Eur. 8.0 ·

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Polysorbate 80	5,0
Lecithin	0,7
Histidine	1,0
Casein Peptone	15,0
Soy Peptone	5,0
Sodium Chloride	5,0
Sodium Thiosulphate	0,5
Agar	15,0
pH: 7.3±0.2	