



**PRODUCT CODE: 456213**

**Sabouraud Glucose Agar+Chloramphenicol (Ph. Eur.) (irradiated) (Prepared Plate (@ 90 mm)) for microbiology**

Physical Data: Prepared media of light beige colour · Procedure: Inoculate the plate with the sample ·

---

TARIC 3821 00 00 00

**SYNONYMS:** Medium C•

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** USP 38 • Ph. Eur. 8.0 • ATLAS, R. M. and L. C. PARKS (1993) , Handbook of Microbiological Media , CRC Press. Inc. London, PASCULA ANDERSON M<sup>a</sup> R<sup>a</sup> (1992) , Microbiología Alimentaria , Díaz de Santos S.A. Madrid, VANDERZANT & SPLITTSTOESSER (1992) , Compendium of Methods for the Microbiological Examination of Food III Ed. , American Public Health Association. Washington D.C. •

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

D(+)-Glucose	40,0
Chloramphenicol	0,05
Casein Peptone	5,0
Meat Peptone	5,0
Agar	15,0
pH: 5.6±0.2	

