



**PRODUCT CODE: 624572**

**Potassium Iodide solution 30% w/v VINIKIT, for wine analysis**

for determination of reducing sugar in wine, according to Rebelein method

---

IK

KI

M.= 166,01

CAS [7681-11-0]

EINECS 231-659-4

TARIC 2827 60 00 00

**PHYSICAL DATA:** liquid, Clear, Colourless, Miscible with water • D 20/4 1,224 •

**BIBLIOGRAPHY:**

**HAZARDOUS:** LD L0 oral rbt 916 mg/kg • LD L0 oral mus 1.862 mg/kg •



H: H372 •

P: P260 • P264 • P270 • P314 • P501 •

**WEIGHT/VOLUME INFORMATION:** 1l~1,224 kg    1kg~0,817 l

**COMPOSITION:**

Potassium Iodide  
Water (s.q.m)

30 g  
100 ml

**SPECIFICATIONS:**

Assay (as KI) (Iodom.) (w/v)  
Identity :  
Identity

29,9 - 30,1 %  
  
passes test

Ed.: 4 . Vig.: 21.06.2017 .

Prod.: 624572